

BETACHEK® Visual

Blood glucose testing without a meter

Clinic use
or Diabetes
Self-testing



Reliable: ✓

Betachek Visual requires no batteries, no cleaning, no control solutions and no coding.

Effective: ✓

Betachek Visual is a proven blood glucose testing solution for effective control of blood glucose.

Affordable: ✓

Betachek Visual test strips are approximately half the cost of meter read test strips¹.

Accurate: ✓

UK government tests show excellent precision and accuracy ($R^2=0.99$) when read in a meter².

ISO15197: ✓

When read visually by trained patients the correlation with lab has been shown to be $R^2=0.973^3$.

If the test strip colour lies between two colours on the scale the result is the average e.g between 110 & 150 = 130

Both units of measure mmol/l and mg/dl

Colour pictures on the label show the three test steps

Each colour block is created using an advanced coating process to improve colour matching

Small steps in the reading scale improves accuracy

Two colour test zone improves reading accuracy

Large measurement range 1.1-44.4 mmol/l (20-800 mg/dl)

Check the colour of an unused strip against the zero to ensure it is in good condition

Available from:

UK NHS Drug Tariff Listed.



Instructional videos - click here



Made in Australia by:

National Diagnostic Products
7-9 Merriwa Street,
Gordon NSW 2072
Australia

Tel: +61 (2) 94184777
Email: info@betachek.com
Web: www.betachek.com

1. UK Drug Tariff 2. UK Medical Device Agency Report 02034. 3. Data on file.

